

Whole mount *in situ* hybridization

Cohn Group, GUDMAP Consortium

Solutions are mixed and prepared in clean, autoclaved (121°C for 45 minutes) glassware or in new RNase/DNase free polypropylene tubes (Corning® 50mL PP Centrifuge Tubes #430921, 2mL Eppendorf® Safe-Lock microcentrifuge tubes #022363352).

Dissection of embryos

1. Prepare 3 tubes for each embryo,
 - i) one with 1ml 4% paraformaldehyde (PFA) / 1X PBS for the embryo,
 - ii) one with 300µl 25mM NaOH to genotype the embryo.
 - iii) one with 500µl DEPC PBS for the extra tissue collected.
2. Dissect the 12.5 dpc embryos in ice-cold 1X PBS.
3. Immediately fix the bottom half of the embryo in 4% PFA / 1X PBS overnight at 4°C.
4. Stick a forelimb to the side of the 25mM NaOH tube. Put remainder of embryo in tube with 500µl DEPC PBS.

Preparation of embryos

1. Wash embryo 2X with 1X PBS for 5 minutes at room temperature on a platform shaker.
2. Dehydrate in a 25%, 50%, 75% & 100% methanol / 1X PBS series. Incubate in each solution for 10 minutes, on ice using a platform shaker.
3. Store embryos in 100% methanol at -20°C until ready to use.

Day 1: Pre-hybridization

Approximately 2ml of solution is used per embryo, per step.

1. Rehydrate embryos in a 75%, 50% & 25% methanol / DEPC treated PBT series. Incubate in each solution for 10 minutes, on ice using a platform shaker.
2. *All Day 1 steps from here until the pre-hybridization step are carried out at room temperature on a platform shaker.*
3. Pre-warm the pre-hybridization solution in a 65°C oven.
4. Wash embryos 2X in DEPC treated PBT for 5 minutes.
5. Bleach in 4% H₂O₂ for 1 hour.
6. Wash 2X in DEPC treated PBT for 5 minutes.
7. Treat embryo with Proteinase K at a concentration of 10µg/ml for 10 minutes.
8. Wash 2X in DEPC treated PBT for 10 minutes.
9. Incubate in postfix for 20 minutes.
10. Wash 2X in DEPC treated PBT for 5 minutes.

11. Incubate the embryos in the pre-warmed pre-hybridization buffer overnight at 65°C in a hybridization oven with shaking.

Day 2: Hybridization buffer with probe

Approximately 2ml of solution is used per embryo, per step.

1. Pre-warm fresh pre-hybridization buffer to 65°C.
2. Thaw probes on ice.
3. Add probe to 2ml hybridization buffer. Probe concentration: 0.1µg/ml
4. Incubate the embryo in the hybridization/probe buffer for 3 days at 65°C in a hybridization oven with shaking.

Day 5: Stringency washes, blocking & anti-digoxigenin immunological staining

Approximately 2ml of solution is used per embryo, per step.

1. Wash 3X in 2X SSC with CHAPS for 30 minutes at 65°C with shaking.
2. Wash 3X in 0.2X SSC with CHAPS for 30 minutes at 65°C with shaking.
3. Wash 2X in KTBT for 10 minutes at room temperature with shaking.
4. Pre-block in KTBT/GS for 3 hours at 4°C on a platform with shaking.
5. At the same time, prepare blocking solution by adding anti-digoxigenin antibody (1:3000 dilution) and incubate at 4°C for 3 hours.
6. Block in pre-absorbed antibody at 4°C overnight.

Day 6: Washes

Approximately 2ml of solution is used per embryo, per step.

1. Wash 4X in KTBT for 1 hour at room temperature on a platform shaker.
2. Wash in KTBT overnight at 4°C on a platform shaker.

Day 7: Colour development

Approximately 2ml of solution is used per embryo, per step.

1. Wash 2X in NTMT at room temperature, shaking, for 15 minutes.
2. Transfer embryos to 24-well plate.
3. Add NBT/BCIP colour (enough to cover the embryo) at room temperature, cover with foil and place on a platform shaker.
4. Check colour at 10 minutes under the microscope. Re-check periodically until colour development is perfect and record the time of development.
5. Wash 6X with KTBT for 15 minutes at room temperature on a platform shaker.
6. Fix in 4% PFA / 1X PBS overnight at 4°C.
7. Transfer to 1X PBS for storage at 4°C.

8. Photograph as soon as possible in petri dishes with 1% agarose as background.

Solutions – make enough for 2ml per embryo for each step in the protocol

1X PBS (1l)

100ml 10X PBS (Fisher Scientific #BP3994)

900ml ddH₂O

Autoclave at 121°C for 45 minutes.

DEPC treated 1X PBS (1l)

100ml 10X PBS (Fisher Scientific #BP3994)

900µl DEPC (Research Products International Corp #D43060-25.0)

Add ddH₂O to make 1l.

Let sit overnight.

Autoclave at 121°C for 45 minutes.

DEPC treated PBT (1l)

1l DEPC treated 1X PBS

1000µl Triton X-100 (MP Biomedicals #194854)

4% PFA / 1X PBS (1l)

100ml 10X PBS (Fisher Scientific #BP3994)

40g PFA (Sigma #P6148)

Add ddH₂O to make 1l.

pH to 7.2.

Vacuum filter.

Store 40ml aliquots at -20°C.

2X SSC (1l)

100ml 20X SSC (Fisher Scientific #BP13254)

900ml ultrapure H₂O

0.2X SSC (1l)

10ml 20X SSC (Fisher Scientific #BP13254)

990ml ultrapure H₂O

10% CHAPS (10ml)

10ml ddH₂O

1g CHAPS (Sigma #C3023)

KTBT (TBST) 500ml

25ml Tris HCl pH 7.5
15ml 5M NaCl
1.5ml 2M KCl
5ml Tween 20 (Fisher Scientific #BP337-100)
Make up to 500ml with ddH₂O

2X SSC with CHAPS (50ml)

50ml 2X SSC
500µl CHAPS (Sigma #C3023)

0.2X SSC with CHAPS (50ml)

50ml 0.2X SSC
500µl CHAPS (Sigma #C3023)

Proteinase K (10µg/ml) 10ml

10ml DEPC treated PBT
10µl Proteinase K at 10mg/ml (Fisher Scientific # BP1700)

Postfix (20ml)

20ml 4% PFA
160µl glutaraldehyde (Fisher Scientific #O2957-1)

Pre-hybridization buffer (1l)

20g blocking powder (Roche #11096176001)
1g CHAPS (Sigma #C3023)
0.5g heparin (Sigma #H3149)
500ml DI formamide (American Bioanalytical #AB00600)
250ml 20X SSC (Fisher Scientific #BP13254)
5ml tRNA (200mg/ml) (Sigma #R6625)
10ml 0.5M EDTA
Make up to 1l with DEPC treated H₂O.
1ml Triton X-100 (MP Biomedicals #194854)
Dissolve overnight at 65°C.
Store in 50ml aliquots at -80°C.

Block (KTBT/GS) 20ml

16ml KTBT
4ml goat serum (Sigma #G6767)

Polyclonal Sheep anti-Digoxin, Digoxigenin antibody (1:3000) 9ml

9ml KTBT/GS
3 µl anti-DIG (Roche #11093274910)

Hybridization buffer with probe

2ml pre-hybridization buffer
1-2 µl probe (Probe concentration: 0.1µg/ml)

NTMT (50ml)

5ml 1M Tris HCl pH 9.5

1ml 5M NaCl
2.5ml 1M MgCl₂
500µl 1% Triton X-100 (MP Biomedicals #194854)
Make up to 50ml with ddH₂O.

Colour (NBT/BCIP) 1ml

2.25µl NBT (75mg/ml in 70% DMF / ddH₂O) (Roche #11585029001)
3.5µl BCIP (50mg/ml in 100% DMF) (Roche #11585002001)
Make up to 1ml with ddH₂O.